INSTRUCTIONS FOR USE

SpO₂ Cables

 TEMPERATURE LIMITS	 CAUTION
 PROTECT FROM MOISTURE	 NON-STERILE
 LOT CODE	 MANUFACTURER
 QUANTITY	 LATEX FREE <small>Not made with natural rubber latex</small>
 EXPIRY	 READ INSTRUCTIONS FOR USE
 MEDICAL DEVICE	 REFERENCE
 SEPARATE COLLECTION AS WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT	 EC ACCORDING TO EU MDR 2017/745

The classification rules may vary from country to country. According to the EU MDR 2017/745 or the **Australian D Regulations**, SpO₂ connection cables are products of **Class I**. According to the rules of **FDA (USA) or Canada**, SpO₂ connection cables are **Class II**.

- The user has to make sure that the SpO₂ sensor connected to the SpO₂ interface cable suits the electromedical device and the type of patients (adult, infant, neonates)
- Please always read the instruction leaflet for the electromedical device before using the SpO₂ interface cable.
- Remove the SpO₂ interface cable to a safe distance from any sources of electromagnetic radiation.
- In the operating room, ensure that all or part of the SpO₂ interface cable is outside of the operating area.
- Ensure that no liquid can reach the connectors' contacts.
- Do not use a cable or part of a cable if there is any risk to the patient (e.g. damaged insulating material).
- The Nissha Medical Technologies SpO₂ interface cable has not been designed to be used in MRI environment.
- If you need to take a cable out of the department for maintenance or verification, this becomes the responsibility of the user department to clean and disinfect the product before shipping or transportation.
- Do not allow any liquid to get inside the various connectors on the interface cable.
- The class and type of protection (BF, CF) against electric shocks are defined by the type of electromedical device to which the SpO₂ interface cable is connected.
- Always read the instruction manual for the particular device and any accessories to be used before placing the device into service.
- Nissha Medical Technologies may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in this instruction manual.

ISO 13485:2016
HS.NisshaMedical.com/IFU

SpO₂ CABLES

Please carefully read the following instructions:

Failure to observe these precautions for use may lead to undesirable medical consequences.

Important note:

This insert is designed to provide guidance for the use and handling of the SpO₂ interface cable which connects the **SpO₂ sensor** used for the **non-invasive** and continuous **measurement of arterial oxygen saturation**, to the electromedical device.

No reference is made to a specific medical technique. The manufacturer declines any responsibility for problems resulting from improper use of the product.

I. IDENTIFICATION / FIELD OF APPLICATION

IDENTIFICATION

The **SpO₂ interface cable**, designed by **Nissha Medical Technologies**, is the component which connects the **SpO₂ sensor** used for the **non-invasive** and continuous **measurement of arterial oxygen saturation**, to the electromedical device.

The **interface cable** can only be used with the electromedical device (pulse oximeter) for which it has been designed. This is shown either on the device or on its packaging. The compatible devices (sensors, fastening equipment, etc.) are listed in the **Nissha Medical Technologies** commercial catalogue (**COMM/DOCU001/018**) which can be viewed and downloaded on the company's website at HS.Nisshamedical.com/ifu.

FIELD OF APPLICATION

The **Nissha Medical Technologies** SpO₂ interface cable can be used with any patient (adult, infant, neonates) wherever the arterial oxygen saturation needs to be measured using the non-invasive method.

 The user has to make sure that the SpO₂ sensor connected to the SpO₂ interface cable suits the electromedical device and the type of patients (adult, infant, neonates)

ANTIBACTERIAL QUALITIES

The envelopes of the medical bioactive **Nissha Medical Technologies** cables carry an additional chemical compost of non toxic material which gives them a slightly acid surface like the human skin pH neighbour 5. Almost all the pathogenic agents responsible for the nosocomial sicknesses are destroyed in this environment.

Nissha Medical Technologies offers a full range of **SpO₂ interface cables**.

If you wish to place any subsequent orders, please use the code number shown either on the **Nissha Medical Technologies** SpO₂ interface cable or on its packaging;

For further information about this product, please either contact Nissha Medical Technologies or visit the company's website at www.NisshaMedical.com.

II. STORAGE / PACKAGING / SYMBOLS

STORAGE

The storage conditions for **SpO₂ interface cables** are as follows:

- Ambient temperature: 10 to 50 °C- 50 to 122 °F

PACKAGING

Nissha Medical Technologies SpO₂ interface cables are packaged individually.

When it is not in use, the **Nissha Medical Technologies SpO₂ interface cable** must be stored in its original packaging in order to prevent any damage which might reduce its service life, performances and/or safety level.

III. PERFORMANCES / RELIABILITY / SAFETY / COMPATIBILITY / MECHANICAL and ELECTRICAL INTEGRITY / SYMBOLS / ALLERGICITY

PERFORMANCES / RELIABILITY

Nissha Medical Technologies SpO₂ interface cables are inspected both during and at the end of the manufacturing process according to protocols drawn up in line with the specifications of current standards and regulations on this type of product.

The final results of the tests, which have been carried out on a group of representative types of cables according to a technical protocol drawn up by **Nissha Medical Technologies**, have been confirmed by an accredited Laboratory.

(Gmed technical report no. 5358010/1T)

They have also undergone clinical testing and assessments.

SAFETY

Nissha Medical Technologies SpO₂ interface cables are designed and manufactured in accordance with the general and special specifications of the relevant current European and/or international standards:

(Standards EN IEC 60601-1 / EN ISO 9919, etc.)

Nissha Medical Technologies SpO₂ interface cables are an integral part of the "applied part" to the patient as defined by the **EN IEC 60601-1** International safety standard.

The **SpO₂ interface cable's safety class, type of protection (BF, CF) and degree of protection** against electric shocks are closely linked to those of the electromedical device to which it is connected.

Low frequency leakage currents, measured in accordance with the recommendations of the current standards which are applicable to this product, are below the authorised values for the highest safety level.

WARNING

- Please always read the instruction leaflet for the electromedical device before using the SpO₂ interface cable.
- Remove the SpO₂ interface cable to a safe distance from any sources of electromagnetic radiation.
- In the operating room, ensure that all or part of the SpO₂ interface cable is outside of the operating area.
- Ensure that no liquid can reach the connectors' contacts.

PLEASE NOTE

Nissha Medical Technologies may not be held liable for any incidents which might occur to the patient, to the user and to any other persons caused by the presence of dangerous electrical currents coming from the electromedical device if there is any fault at all.

COMPATIBILITY

On the company's website (www.NisshaMedical.com), Nissha Medical Technologies offers its customers a downloadable document including information on the device's compatibility along with technical details about it.

2/48

MECHANICAL and ELECTRICAL INTEGRITY

In order to ensure that the SpO₂ interface cable offers good mechanical resistance to traction and to flexion and to reduce any leakage currents which might travel through the patient, **Nissha Medical Technologies** uses high quality materials for the manufacturing of the devices.

The moulded connectors are fitted with soft inserts designed to minimise the risk of the cable breaking at this point.

The contact surface of the connector pins is treated in such a way as to reduce the contact resistance between two pins as far as possible even after a large number of connections and disconnections.

(Test report no. CPB1097 – CPB1897/CPB2897)

Nissha Medical Technologies SpO₂ interface cables can withstand repeated defibrillation shocks. They have no accessible metal parts.

When it is not in use, a cable must be stored in its original packaging in order to prevent any untimely damage to its insulating sheath and its conductors which might reduce its lifetime, performances and/or safety level.

SAFETY SYMBOLS

Explanation of the symbols used on the electromedical device:

or

Shows that the SpO₂ cables and accessories are fitted with special protection against electric shocks (including admissible leakage currents) and defibrillation shocks.

The class and type of protection (BF, CF) against electrical shocks are defined by the type of electromedical device to which the Nissha Medical Technologies ECG cable is connected to.

ALLERGICITY

The insulating materials used to manufacture **Nissha Medical Technologies SpO₂ interface cables** have been subjected to allergicity tests. These tests have not shown the presence of any products which might develop an intolerable allergic reaction.

IV. INSTALLATION / USE / MAINTENANCE / HYGIENE / STERILISATION

INSTALLATION


Please follow these instructions for optimum installation and operation of the **SpO₂ interface cable**:

(see also the instruction manual for the particular electromedical device)

On the patient's side:

- First of all connect the interface cable to the appropriate **SpO₂ sensor**.
- Place the **SpO₂ sensor** on the chosen site on the patient, which must be suitable for the measurement, referring to the instructions in the leaflet which came with the particular electromedical device. Make sure the SpO₂ sensor suits the patient (adult, infant, neonate)

Example below showing how to position the sensor on the finger:


Do not push the finger in too deeply

On the device's side:

- Connect the **SpO₂ interface cable** to the appropriate device (pulse oximeter).
- Start the device and make sure that it is working properly by reading the instruction leaflet for the electromedical device.

Examples of plugs:

3/48

HOW TO USE:

GENERAL CONDITIONS

The conditions for using the **Nissha Medical Technologies SpO₂ interface cable** are as follows:

- Ambient temperature: 10 to 50 °C- 50 to 122 °F

SPECIAL CONDITIONS

- **Do not use an SpO₂ interface cable if there is any risk to the patient (e.g. damaged insulating material).**

PREVENTIVE MAINTENANCE

The service life of the **SpO₂ interface cable** depends on a large number of parameters e.g.

- How many times it is applied
- Adherence to preventive maintenance
- Being kept in a good state of hygiene

A regular visual and electrical inspection of the conductors will determine whether the **SpO₂ interface cable** needs to be replaced.

WARNING:

 **If you need to take a cable out of the department to another department for maintenance or expertise, it is the responsibility of the department which has made use of the cable to clean and disinfect it before transferring or moving it.**

CORRECTIVE MAINTENANCE

- There is no corrective maintenance for this product. If it is damaged it must be discarded immediately.

HYGIENE

PLEASE NOTE:

Do not allow any liquid to get inside the various connectors on the interface cable.

When using the **Nissha Medical Technologies SpO₂ interface cable** in departments where a level of sterility must be adhered to, you must proceed to use the cleaning and/or disinfection methods described below.

RE-PROCESSING INSTRUCTIONS OF SpO₂ INTERFACE CABLES (Acc. to EN ISO 17664:2004)

Manufactured by: Nissha Medical Technologies	
Devices: SpO ₂ Cables	
CAUTION	When cleaning or disinfecting the cables make sure that the plugs at the end of the cables do not become immersed in any liquid to avoid any electrical problems
Re-processing limits	Nissha Medical Technologies SpO ₂ interface cables are designed to support repetitive defibrillation shocks and cleaning and disinfection cycles. Any misuse of the cable may shorten its lifetime.
INSTRUCTIONS	
In the place of use:	Remove any stains or dirt with a soft disposable cloth
Isolation and transportation:	Remember to always reprocess the device after use
Cleaning preparation:	No specific requirement

Automatic cleaning:	Is strictly not applicable
Manual cleaning:	<ul style="list-style-type: none"> • Clean the Nissha Medical Technologies SpO₂ interface cable and its connecting cable (including the connector) with a cloth soaked in soapy warm water
Disinfection:	<p style="text-align: right;">(Taken from the study by the ANIOS Laboratory, no. 6416.94/0387)</p> <p>Method A</p> <ul style="list-style-type: none"> • fill a dip tank with a 0.5% HEXANIOS G+R solution • partially immerse the cable keeping the ends dry in order to avoid any electrical problems on the connectors • allow a contact time of 15 minutes • rinse the cable • dry the cable with absorbent paper • during the immersion rub the ends of the cable using a LINGET ANIOS wipe. <p style="text-align: right;">(Replace the solution every 48 hours) (Taken from the study by the ANIOS Laboratory, no. 14496.02/052)</p> <p>Method B (new cold method)</p> <ul style="list-style-type: none"> • fill a dip tank with an ANIOXYDE 1000 activated preparation • but do not use before 30 minutes • partially immerse the cable protecting the ends in order to avoid any electrical problems on the connectors • allow a contact time of 10 to 30 minutes depending upon the required level of disinfection • regularly check the peracetic acid content using a test strip • and afterwards rinse with ANIOS mains water (pH = 7.3; TH = 48°F)
Drying:	Dry the cable with absorbent paper
Maintenance, control and tests:	Please check visually that the cable is in working order
Conditionment	Please follow the standard protocol of packaging of your institution before sterilization
Sterilization:	<ul style="list-style-type: none"> • sterilization time: 22 hours which takes into account 20 hours of gas exposure • sterilizing agent: Carbon Dioxide + Ethylene Oxide (80/20%) • initial vacuum: -70 kPa • relative humidity: >60% • temperature: 50 °C • Ethylene oxide concentration: 530 g/m³ • pressure: 120 kPa • final vacuum: -70 kPa • aeration <p>(Test report 001: sterilization test / ECG cables – RE/IP/BPF – Revision 0 dated 26/02/98)</p>
Please Note:	Never sterilize the cables and accessories in an "autoclave" using boiling water or steam
Preservation:	Please read the conditions of storage in this document
Complementary information:	The above methods have been validated and are strongly advisable. Particularly we advise you strongly to use the recommended or equivalent products above and not to exceed the indicated time of application or immersion; otherwise the life time of the cable might be reduced.
Contact the manufacturer:	See below or visit the website www.NisshaMedical.com

V. GUARANTEE/LIABILITY

Any unused **SpO₂ interface cable** kept in its original packaging and which has not suffered any visible damage is guaranteed by **Nissha Medical Technologies** for one year.

Nissha Medical Technologies guarantees the device's conformity to the specifications of the current safety and performance standards applicable to it.

PLEASE NOTE:

The class and type of protection against electric shocks are also linked to those of the electromedical device (pulse oximeter) to which the sensor and its cable are connected.

Always read the instruction manual for the particular device and any accessories used with it before placing the device into service.

Nissha Medical Technologies may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in this instruction manual.

NISSHA
MEDICAL TECHNOLOGIES

IP INTEGRAL
PROCESS
Cables & Leadwires

MODE D'EMPLOI

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 TEMPERATURE LIMITS	 CAUTION
 PROTECT FROM MOISTURE	 NON-STERILE
 LOT CODE	 MANUFACTURER
 QUANTITY	 LATEX FREE
 EXPIRY	 READ INSTRUCTIONS FOR USE
 MEDICAL DEVICE	 REFERENCE
 SEPARATE COLLECTION AS WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT	 EC ACCORDING TO EU MDR 2017/745

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RE-PROCESSING INSTRUCTIONS OF SpO₂ INTERFACE CABLES (Acc. to EN ISO 17664:2004)

Manufactured by: Nissha Medical Technologies	
Devices: SpO ₂ Cables	
CAUTION	No dejar que penetre líquido en el interior de los conectores del cable de interfaz.
Re-processing limits	Nissha Medical Technologies SpO ₂ interface cables are designed to support repetitive defibrillation shocks and cleaning and disinfection cycles. Any misuse of the cable may shorten its lifetime.
INSTRUCTIONS	
In the place of use:	Remove any stains or dirt with a soft disposable cloth
Isolation and transportation:	Remember to always reprocess the device after use
Cleaning preparation:	No specific requirement

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
NISSHA
MEDICAL TECHNOLOGIESIP INTEGRAL
PROCESS
Cables & Leadwires

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Important note:

This insert is designed to provide guidance for the use and handling of the SpO₂ interface cable which connects the SpO₂ sensor used for the **non-invasive** and continuous **measurement of arterial oxygen saturation**, to the electromedical device.

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TEMPERATURE LIMITS	CAUTION
PROTECT FROM MOISTURE	NON-STERILE
LOT LOT CODE	MANUFACTURER
QTY QUANTITY	LATEX FREE
EXPIRY	READ INSTRUCTIONS FOR USE
MEDICAL DEVICE	REFERENCE
SEPARATE COLLECTION AS WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT	EC ACCORDING TO EU MDR 2017/745

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The user has to make sure that the SpO₂ sensor connected to the SpO₂ interface cable suits the electromedical device and the type of patients (adult, infant, neonates)

- Please always read the instruction leaflet for the electromedical device before using the SpO₂ interface cable.
- Remove the SpO₂ interface cable to a safe distance from any sources of electromagnetic radiation.
- In the operating room, ensure that all or part of the SpO₂ interface cable is outside of the operating area.
- Ensure that no liquid can reach the connectors' contacts.
- Do not use a cable or part of a cable if there is any risk to the patient (e.g. damaged insulating material).
- The Nissha Medical Technologies SpO₂ interface cable has not been designed to be used in MRI environment.
- If you need to take a cable out of the department for maintenance or verification, this becomes the responsibility of the user department to clean and disinfect the product before shipping or transportation.
- Do not allow any liquid to get inside the various connectors on the interface cable.
- The class and type of protection (BF, CF) against electric shocks are defined by the type of electromedical device to which the SpO₂ interface cable is connected.
- Always read the instruction manual for the particular device and any accessories to be used before placing the device into service.
- Nissha Medical Technologies may not be held liable for any incidents which might occur in the event of any failure to adhere to the rules of installation and use mentioned in this instruction manual.

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Important note:

This insert is designed to provide guidance for the use and handling of the SpO₂ interface cable which connects the SpO₂ sensor used for the non-invasive and continuous measurement of arterial oxygen saturation, to the electromedical device.

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The user has to make sure that the SpO₂ sensor connected to the SpO₂ interface cable suits the electromedical device and the type of patients (adult, infant, neonates)

ANTIBACTERIAL QUALITIES

The envelopes of the medical bioactive Nissha Medical Technologies cables carry an additional chemical compost of non toxic material which gives them a slightly acid surface like the human skin pH neighbour 5. Almost all the pathogenic agents responsible for the nosocomial sicknesses are destroyed in this environment.

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Complementary information:	The above methods have been validated and are strongly advisable. Particularly we advise you strongly to use the recommended or equivalent products above and not to exceed the indicated time of application or immersion; otherwise the life time of the cable might be reduced.
Contact the manufacturer:	See below or visit the website www.NisshaMedical.com

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TEMPERATURE LIMITS	CAUTION
PROTECT FROM MOISTURE	NON-STERILE
LOT CODE	MANUFACTURER
QUANTITY	LATEX FREE
EXPIRY	READ INSTRUCTIONS FOR USE
MEDICAL DEVICE	REFERENCE
SEPARATE COLLECTION AS WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT	EC ACCORDING TO EU MDR 2017/745

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RE-PROCESSING INSTRUCTIONS OF SpO₂ INTERFACE CABLES (Acc. to EN ISO 17664:2004)

Manufactured by: Nissha Medical Technologies	
Devices: SpO ₂ Cables	
CAUTION	Ervoor zorgen dat er geen vloeistoffen binnendringen in de verschillende schakelaars van de verbindingkabel.
Re-processing limits	Nissha Medical Technologies SpO ₂ interface cables are designed to support repetitive defibrillation shocks and cleaning and disinfection cycles. Any misuse of the cable may shorten its lifetime.
INSTRUCTIONS	
In the place of use:	Remove any stains or dirt with a soft disposable cloth
Isolation and transportation:	Remember to always reprocess the device after use
Cleaning preparation:	No specific requirement

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RE-PROCESSING INSTRUCTIONS OF SpO₂ INTERFACE CABLES (Acc. to EN ISO 17664:2004)

Manufactured by: Nissha Medical Technologies	
Devices: cabo de interface SpO ₂	
Atenção	Não fazer com que o líquido penetre no interior dos conectores do cabo de interface.
Re-processing limits	Nissha Medical Technologies SpO ₂ interface cables are designed to support repetitive defibrillation shocks and cleaning and disinfection cycles. Any misuse of the cable may shorten its lifetime.
INSTRUCTIONS	
In the place of use:	Remove any stains or dirt with a soft disposable cloth
Isolation and transportation:	Remember to always reprocess the device after use
Cleaning preparation:	No specific requirement

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