

End Tidal Sampling: EtCO₂ /O₂ Cannulas

FEATURES



RELIABLE ACCURATE VALUES

Provides clean, crisp and real-time waveforms comparable to leading capnography sampling lines



DIVIDED CANNULA

Enables reliable O₂ delivery and CO₂ sampling



PATIENT COMFORT

Soft, pliable headset tubing for exceptional comfort



MULTIPLE SIZES AND APPLICATIONS


Available in adult and pediatric sizes



ECONOMICAL SOLUTION

High value performance compared to leading ETCO₂ sampling lines

Non-invasive real-time breath-by-breath examination of a patient's ventilatory condition



700000360

Adlt. Oral/Nasal
Divided Cannula, 7' O₂, 6' CO₂ line



Reflective Connector for Extended Use



Reflective Connector for Short-Term Use




End Tidal Sampling: EtCO₂ /O₂ Cannulas

SALTER LABS PRODUCT LINE - 25/case






Inline Sampling Lines with 15mm Adapter



Part #	Reference	Description	Size	Interface	O ₂ Length	CO ₂ Length	
700000356	4MSF5-6-25	EtCO ₂ Filtered Sampling Line with Adapter	Adult	-	-	6'	
700000354	4MSF5-13-25	EtCO ₂ Filtered Sampling Line with Adapter	Adult	-	-	13'	
700000355	4MSF5-L-6-25	EtCO ₂ Filtered Sampling Line with Adapter, Nafion™ Tubing	Adult	-	-	6'	





Nasal CO₂ Sampling with O₂ Delivery



700000357	4MSF3-7-6-25	EtCO ₂ /O ₂ Filtered Capnography Cannula	Adult	Nasal	7'	6'	
700000359	4MSF3-L-7-6-25	**MOST POPULAR EtCO ₂ /O ₂ Filtered Capnography Cannula, Nafion™ Tubing	Adult	Nasal	7'	6'	
700000362	4MSF3-L-PED-7-6-25	EtCO ₂ /O ₂ Filtered Capnography Cannula, Nafion™ Tubing	Pediatric	Nasal	7'	6'	

Oral/Nasal CO₂ Sampling with O₂ Delivery



700000358	4MSF1-7-6-25	EtCO ₂ /O ₂ Filtered Capnography Cannula	Adult	Oral/Nasal	7'	6'	
700000360	4MSF1-L-7-6-25	**MOST POPULAR EtCO ₂ /O ₂ Filtered Capnography Cannula, Nafion™ Tubing	Adult	Oral/Nasal	7'	6'	
700000361	4MSF2-PED-7-6-25	EtCO ₂ /O ₂ Filtered Divided Cannula	Pediatric	Oral/Nasal	7'	6'	
700000363	4MSF2-L-PED-7-6-25	EtCO ₂ /O ₂ Filtered Divided Cannula, Nafion™ Tubing	Pediatric	Oral/Nasal	7'	6'	

*The OEM trademarks identified herein are the trademarks of the respective OEMs, and not of Nissha Medical Technologies. Nissha Medical Technologies disclaims any affiliation, connection, or association between its products and those of the respective OEMs.