

Medical Products

Today's healthcare providers are faced with a variety of demands and challenges. Nissha Medical Technologies offers solutions geared toward the needs of both patient and provider with fast, efficient and easy-to-use products with a priority on function and patient comfort.

ECG Electrodes | Defibrillation Pads | Patient Monitoring Kits
Patient Monitoring Charts | Ultrasound Film | Grounding Pads



ECG ELECTRODES




In critical diagnostic situations, use Nissha Medical Technologies ECG electrodes to ensure quality readings, zero artifact and no repeat procedures.




Defibrillation • Cardioversion • Pacing • ECG Monitoring

DEFIBRILLATION PADS



HeartSync Defibrillation Electrodes are designed with your needs in mind — featuring long leadwires, 36 month shelf life, and ease of placement.



CLARAVUE® PATIENT MONITORING KITS



CLARAVUE® is the best-in-class, hygienic solution for simple, accurate and fully radiolucent patient monitoring.



DISPOSABLE,
HYGIENIC
PREWIRED
ELECTRODES



QUICK
PATIENT
HOOKUP



RADIOLUCENT
FOR
UNOBSTRUCTED
IMAGING



PRE-
POSITIONED
ELECTRODE
PAIRS

NISSHA MEDICAL TECHNOLOGIES ULTRASOUND FILM



NMT is value leader in ultrasound film. High-density and high-gloss formats available.




HIGH
QUALITY
IMAGING



WATER
RESISTANT

SWAROPLATE® SPLIT | NON-SPLIT | TWIN-SAFE | OVALPLATE

GROUNDING PADS



From our patented TwinSafe design to our standard SwaroPlate offering, Grounding (Diathermy) Pads are compatible with many electrosurgical generators.




COMPLETE
CONTACT
DURING
HIGH
FREQUENCY
SURGERY



ROUNDED
CONDUCTIVE
PLATE FOR
REDUCED
EDGE EFFECT
BURN RISK

Graphic Controls by NISSHA

CHART PAPER



NMT is the ORIGINAL Chart Paper manufacturer, offering the widest assortment of chart papers for all departments in the hospital.



PRECISE
GRID
MARKINGS



MEETS
OEM
SPECS