



# Cardiology & Stress Lab

CLARAVUE® provides the best-in-class, hygienic solution for simple, accurate, and fully radiolucent patient monitoring. CLARAVUE® includes disposable pre-wired and pre-positioned electrodes and patented reusable bioactive trunk cables.

For the past ten years, CLARAVUE® has performed as the proven choice for superior trace quality and infection prevention in prominent European hospitals.

Lightweight leadwires (with a direct connection to the electrode sensor) help to reduce artifact

## FEATURES



### PRE-WIRED ELECTRODES

Help to reduce motion artifact.



### TWO SIMPLE & QUICK CONNECTIONS

Consolidate 10 leadwire ends into 2 simple connections.



### PRE-POSITIONED ELECTRODE PAIRS

Increase efficiency, saving time & money.



### PATENTED CONNECTOR CABLE

Standardized for all monitors.




### FULLY RADIOLUCENT LEADWIRES

Reduce clutter in imaging.



### HYGIENIC, DISPOSABLE SOLUTION:

Reduces risk of cross contamination.



## PACKAGING

Electrodes are packaged in a single-use pouch with an easy-to-open tear.

### Box




### Single-Use Kit




FRONT

BACK






50601-US



50602-US



50604-US

## CLARAVUE® RADIOLUCENT DISPOSABLE PRE-WIRED ELECTRODES

Ref	Product #	Description	Per Kit	Kits/Box	Kits/Case	Solid Gel	Duration of Use
50601-US	32028778	Adult Patient Kit - 5 electrodes	5	30	360		
50602-US	32028779	Adult Patient Kit - 10 electrodes	10	25	300		
50604-US	32028781	Adult Patient Kit - 5 electrodes - V-Series	5	30	360		

## COMPATIBLE MONITOR CONNECTIONS

Compatible cables include adaptors for Bard, Biosense, Mindray, Datascope, GE Marquette, Philips, Datex, Mennen, Mortara, Nihon Kohden, GE, Getemed, Spacelabs, Siemens, Drager, Schiller, Zoll, and more.



Mindray



Philips



Spacelabs  
TruLink



GE



Physio



Nihon  
Kohden

The OEM trademarks identified herein are the trademarks of the respective OEMs, and not Vermed. Vermed disclaims any affiliation, connection, or association between its products and those of the respective OEMs.

## CABLE CONNECTIONS

ELECTRODE → PATENTED ELECTRODE CABLE CONNECTOR → CABLE → MONITOR

