



## Defibrillator Compatibility

Match Defibrillator Brands to Vermed Part Numbers



BRAND OF DEFIBRILLATOR	MODEL OF DEFIBRILLATOR	VERMED COMPATIBLE DEFIBRILLATION ELECTRODE PART NUMBERS	
<b>Philips Medical</b>	<ul style="list-style-type: none"> <li>• XLT</li> <li>• HeartStart XL</li> <li>• HP Codemaster XL (Plug Style Connector)</li> </ul>	<ul style="list-style-type: none"> <li>• XL+</li> <li>• MRX</li> <li>• MRX ALS</li> </ul>	<ul style="list-style-type: none"> <li>• D311</li> <li>• D309</li> <li>• D331</li> </ul>
	<ul style="list-style-type: none"> <li>• HeartStart FR2+</li> </ul>	<ul style="list-style-type: none"> <li>• HeartStart Forerunner</li> </ul>	<ul style="list-style-type: none"> <li>• D311</li> <li>• D309</li> </ul>
	<ul style="list-style-type: none"> <li>• HeartStart FR3</li> </ul>	<ul style="list-style-type: none"> <li>• HeartStart FRx</li> </ul>	<b>NOT COMPATIBLE</b>
<b>Physio-Control</b>	<ul style="list-style-type: none"> <li>• LP-9P</li> <li>• LP-10C</li> <li>• LP-12</li> </ul>	<ul style="list-style-type: none"> <li>• LP-15</li> <li>• LP-20</li> </ul>	<ul style="list-style-type: none"> <li>• D304</li> <li>• D313</li> <li>• D319</li> </ul>
	<ul style="list-style-type: none"> <li>• LP-500</li> <li>• LP-1000</li> </ul>	<ul style="list-style-type: none"> <li>• LP-500 DPS</li> </ul>	<ul style="list-style-type: none"> <li>• D304</li> <li>• D313</li> </ul>
<b>Zoll Medical</b>	<ul style="list-style-type: none"> <li>• M-Series</li> <li>• R-Series</li> </ul>	<ul style="list-style-type: none"> <li>• E-Series</li> <li>• X-Series</li> </ul>	<ul style="list-style-type: none"> <li>• D303</li> <li>• D312</li> <li>• D314</li> </ul>
	<ul style="list-style-type: none"> <li>• Zoll AED PLUS</li> </ul>		<b>NOT COMPATIBLE</b>
<b>Cardiac Science</b>	<ul style="list-style-type: none"> <li>• PowerHeart AED G3 Plus</li> <li>• PowerHeart AED G3 Pro</li> </ul>		<ul style="list-style-type: none"> <li>• D326</li> </ul>

Please contact Vermed directly for any questions regarding compatibility.